

Tube: Stationary Anode Tubes – Medical, Dental, & Veterinary Use

Model Numbers and Specifications

Tube Type	Max Operating Voltage (KV)	Nominal Focal Spot Value	Auod Heat Storage Capacity (KG)	Max Rating (KW) (0.1S)	Filament Max Current (A)	Data Max Voltage (V)	Target Angle (degree)	Outline Diameter (mm)	Dimension Length (mm)	Remarks
XD1-3/100	100	2.3	30	3	4.5	7	19	50	165	
XD2-1.4/85	85	1.5	9	1.4	2.5	3.3	23	35	110	
XD3-3.5/100	100	2.6	35	3.5	4.5	7.8	19	50	180	
XD4-2.9/100	100 100	1.80 4.3	60 60	2 9	4.5 4.5	5.5 11	19 19	60 60	247 247	
XD6-1.1, 3.5/100	100 100	1.0 2.6	35 35	1.1 3.5	4.5 4.5	4.3 7.8	19 19	50 50	180 180	
XD7-1.05/35	35	1	35	1.05	4.0	4.2	19	55	148	Target Material Mo
XD8-0.1/70	70	0.3	4	0.1	1.5	3.0	12	26	80	
XD9-4/100	100	2.2	40	4	4.0	7.5	19	50	147	
XZ1-4/250	250	6	-	4	4.0	6.6	30	105	446	Deep Part Therapy Tube with Grid Rg 30k-50k
X763*	70	0.7	7	0.84	3.3	4.0	12	31	80	
X714A	100 100	0.8 4.2	35 35	1.0 8	4.5 4.5	4.3 11	19 19	50 50	180 180	
X714B	120 120	1.0 2.6	32 32	1.1 3.5	4.5 4.5	4.3 7.8	19 19	60 60	185 185	
X714C	110 100	0.8 2.6	40 40	0.8 3.6	4.5 4.5	4.3 7.8	19 19	60 60	205 205	
X762	180	0.8	30	-	4.4	5.3	32	68	188	Target Material W
X777	90	1.8	35	0.25	4.0	6.5	19	46	145	Replace DE 90
X780 Panorama*	90	0.5	30	1.53	4.0	7.5 2.75	5	45	152	Replace SR90/15F
X782	110 110	0.6 1.5	30 30	0.5 3.2	4.5 4.5	5 7	- -	50 50	165 165	
X782D	110 110	0.6 1.5	30 30	0.5 3.2	4.5 4.5	5 7	15 15	50 50	147 147	
X784*	70	0.25	7	0.2	2.5	2.2	12	30	75	
S1	60	0.4	10	0.1	2.2	3.3	12	40	110	
XD12*	70	0.8	7	0.56	2	3.5	19	30.5	73	
X790	110 110	1 4.2	30 30	0.8 3.6	4.5 4.5	4.2 11	19 19	51 51	197 197	
X791	3.5 3.5	0.5 1	- -	- 1.05	4 4	4.3 4.2	19 19	55 55	148 148	

*Dental Tube

Tube: Stationary Anode X-Ray Tube for Treatment

Tube Type	Tube Voltage (kV)	Nominal Focal Spot	Angle Degrees	Nominal Anode Input Power (kW)	Filament Data		Data Max Voltage (V)	
					Max. Current (A)	Max Voltage (V)	Diameter	Length
XD1-4/250	250	6	30	4 (15 min)	4	6 +/- 0.6	105	445